



Announcing the World's Smallest High-Resolution Manometry Catheter...

Provides Maximum Patient Comfort and Diagnostic Accuracy

- Advanced Sensing Technology**
 The small diameter High-Resolution catheter incorporates the latest in sensing technology for simplified and reliable data acquisition
- Complete Coverage of Entire Esophagus**
 36 pressure sensors provide coverage from the pharynx to the stomach with single placement of the catheter
- True Circumferential Sensing**
 Each channel contains 12 individual pressure data points
- Improved Patient Comfort**
 Small 2.7 diameter (7.5 French) improves ease of intubation and increases patient tolerance



Pediatric Small Diameter Catheter

- 2.7 mm in diameter (7.5 French) spaced 7.5 mm apart
- Manoshield™ disposable catheter sheath provides sanitary protection

Adult Small Diameter Catheter

- 2.7 mm in diameter (7.5 French) spaced 1 cm apart
- Manoshield™ disposable catheter sheath provides sanitary protection

